



Pediatric Surgical Associates, Ltd.

Patient and Family Education

ACE Stopper

Medicina's ACE stoppers are made of implantable grade silicone and are used to plug/seal stoma sites at various parts of the body. ACE stoppers are soft, comfortable and come in a variety of lengths and widths depending on the size of the stoma. Stoppers should be discarded and replaced every month or sooner if there is any damage to the joint between the skin disc and the stopper shaft. They all have a circular 15mm skin disc and can perform two functions:


- ACE stoppers form a seal in between insertions of the catheter if the stoma is leaking
- Ace stoppers maintain the patency of the stoma if they catheter insertion is difficult or the site closes in

ACE Dressing 30 per box

8-14FR diameter Keyhole Dressing 30 per box

16-20FR diameter Keyhole Dressing 30 per box

*Latex free, Non-DEHP



8FR Ace Stopper, Length 15mm

8FR Ace Stopper, Length 30mm

8FR Ace Stopper, Length 60mm

10FR Ace Stopper, Length 15mm

10FR Ace Stopper, Length 30mm

10FR Ace Stopper, Length 60mm

12FR Ace Stopper, Length 30mm

12FR Ace Stopper, Length 60mm

12FR Ace Stopper, Length 100mm

14FR Ace Stopper, Length 30mm

14FR Ace Stopper, Length 60mm

14FR Ace Stopper, Length 100mm

*Latex free, Non-DEHP