



Pediatric Surgical  
Associates, Ltd.

Patient and Family Education


## ACE Stopper

Medicina's ACE stoppers are made of implantable grade silicone and are used to plug/seal stoma sites at various parts of the body. ACE stoppers are soft, comfortable and come in a variety of lengths and widths depending on the size of the stoma. Stoppers should be discarded and replaced every month or sooner if there is any damage to the joint between the skin disc and the stopper shaft. They all have a circular 15mm skin disc and can perform two functions:

- ACE stoppers form a seal in between insertions of the catheter if the stoma is leaking
- Ace stoppers maintain the patency of the stoma if they catheter insertion is difficult or the site closes in

ACE Dressing	30 per box
8-14FR diameter Keyhole Dressing	30 per box
16-20FR diameter Keyhole Dressing	30 per box

\*Latex free, Non-DEHP



8FR Ace Stopper, Length 15mm
8FR Ace Stopper, Length 30mm
8FR Ace Stopper, Length 60mm
10FR Ace Stopper, Length 15mm
10FR Ace Stopper, Length 30mm
10FR Ace Stopper, Length 60mm
12FR Ace Stopper, Length 30mm
12FR Ace Stopper, Length 60mm
12FR Ace Stopper, Length 100mm
14FR Ace Stopper, Length 30mm
14FR Ace Stopper, Length 60mm
14FR Ace Stopper, Length 100mm

\*Latex free, Non-DEHP